

Product datasheet for TA336770

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

beta Actin (ACTB) Mouse Monoclonal Antibody [Clone ID: 8H10D10]

Product data:

Product Type: Primary Antibodies

Clone Name: 8H10D10

Applications: ELISA, FC, ICC/IF, Immunoblotting, WB

Recommended Dilution: Flow (Intracellular): 1:200 - 1:400, Western Blot: 1:1000-1:5000, Simple Western: 12.5, Flow

Cytometry: 1:200-1:400, ELISA: 1:10000, Immunocytochemistry/ Immunofluorescence: 1:200-

1:1000, Immunoblotting

Reactivity: Human, Mouse, Rat, Primate, Hamster

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Synthetic peptide corresponding to amino-terminal residues of human Beta Actin,

conjugated to KLH. [UniProt# P60709]

Formulation: PBS, 0.03% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid

freeze-thaw cycles.

Concentration: lot specific

Purification: Ammonium sulfate precipitation

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 42 kDa

Gene Name: actin, beta

Database Link: NP 001092

Entrez Gene 11461 MouseEntrez Gene 81822 RatEntrez Gene 60 Human

P60709



beta Actin (ACTB) Mouse Monoclonal Antibody [Clone ID: 8H10D10] - TA336770

Background: Beta Actin is one of six different actin isoforms which have been identified. Actins are highly

conserved proteins that are involved in cell motility, structure, and integrity. Because Beta Actin is ubiquitously expressed in all eukaryotic cells, it is frequently used as a loading control for assays involving protein detection, such as Western blotting. Antibodies to Beta Actin provide a specific and useful tool in studying the intracellular distribution of Beta Actin and

the static and dynamic aspects of the cytoskeleton.

Synonyms: BRWS1; PS1TP5BP1

Note: This Beta Actin (8H10D10) antibody is useful in Western blot, Flow Cytometry,

Immunocytochemistry/Immunofluorescence and ELISA.


Protein Families: ES Cell Differentiation/IPS

Protein Pathways: Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated

cardiomyopathy, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Leukocyte transendothelial migration, Pathogenic Escherichia coli infection, Regulation of actin


cytoskeleton, Tight junction, Vibrio cholerae infection, Viral myocarditis

Product images:




Simple Western: beta-Actin Antibody (8H10D10) TA336770 - Image shows a specific band for Beta Actin in 1.0 mg/mL of HeLa lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.







Immunocytochemistry/Immunofluorescence: beta-Actin Antibody (8H10D10) TA336770 - Beta actin was detected in NIH-3T3 cells fixed with methanol using mouse anti-mouse beta-Actin monoclonal antibody (TA336770) at 1:5400 dilution. Cells were stained using NorthernLights 557 conjugated anti-mouse secondary antibody (NL007) and counterstained with DAPI.




Flow Cytometry: beta-Actin Antibody (8H10D10) TA336770 - Analysis of MCF-7 cells using anti-Beta Actin mAb (right) and negative control (left).




Immunocytochemistry/Immunofluorescence: beta-Actin Antibody (8H10D10) TA336770 - HeLa cells were fixed and permeabilized for 10 minutes with -20C MeOH. The cells were incubated with beta-Actin Antibody [8H10D10] conjugated to Biotin (TA336770B) at 5ug/ml overnight at 4C and detected with Streptavidin conjugated to DyLight 550. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 100X objective and digitally deconvolved.



Flow Cytometry: beta-Actin Antibody (8H10D10) TA336770 - Analysis of HeLa cells using mouse Monoclonal beta-Actin antibody (Orange) and Isotype control Antibody (Blue).





Flow Cytometry: beta-Actin Antibody (8H10D10) - BSA Free TA336770 - An intracellular stain was performed on NIH3T3 cells with beta-Actin Antibody (8H10D10) TA336770 (blue) and a matched isotype control MAB004 (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 1 ug/mL for 30 minutes at room temperature, followed by Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody, Dylight 550 (84540, Thermo Fisher).




Flow Cytometry: beta-Actin Antibody (8H10D10) - BSA Free TA336770 - An intracellular stain was performed on U-251 MG cells with beta-Actin Antibody (8H10D10) TA336770 (blue) and a matched isotype control MAB004 (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 1 ug/mL for 30 minutes at room temperature, followed by Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody, Dylight 550 (84540, Thermo Fisher).










Western Blot: beta-Actin Antibody (8H10D10) TA336770 - Analysis using anti-Beta Actin mAb against human, mouse, hamster, rat, and primate cell lystate(s): (1) NIH/3T3, (2) Jurkat, (3) HeLa, (4) CHO, (5) PC12, (6) HEK293, (7) COS, (8) A549, and (9) MCF-7. A specific band can be observed at a molecular weight of approximately 43 kDa in all lanes.



Western Blot: beta-Actin Antibody (8H10D10) TA336770 - Lysates of HeLa human cervical epithelial carcinoma cell line, A431 human epithelial carcinoma cell line, K562 human chronic myelogenous leukemia cell line, NIH-3T3 mouse embryonic fibroblast cell line, and C6 rat glioma cell line were probed with 1:2500 mouse anti-beta-Actin monoclonal (TA336770) followed by 1:2000 dilution of donkey anti-mouse IgG-HRP secondary antibody (HAF018).