

Product datasheet for TA336908

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

TrkC (NTRK3) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ICC/IF, IHC

Recommended Dilution: Immunohistochemistry: 1:200, Immunocytochemistry/ Immunofluorescence: 1:1000,

Immunohistochemistry-Paraffin: 1:200

Reactivity: Human, Mouse

Host: Rabbit

Clonality: Polyclonal

Immunogen: A synthetic peptide made to the extracellular domain of the human TrkC protein (within

residues 300-400). [UniProt Q16288].

Formulation: PBS, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid

freeze-thaw cycles.

Concentration: lot specific

Purification: Immunogen affinity purified

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: neurotrophic receptor tyrosine kinase 3

Database Link: NP 002521

Entrez Gene 18213 MouseEntrez Gene 4916 Human

Q16288





Background:

Neurotrophic tyrosine kinase receptor type 3 (NTRK3), commonly known as TRKC, is a member of the protein kinase superfamily, the tyrosine (Tyr) protein kinase family, and more specifically, the insulin receptor subfamily. TRKC is a single pass type I membrane protein that phosphorylates itself (via ligand-mediated auto-phosphorylation) and other proteins in the MAPK pathway upon neurotrophin binding. Mutations to the TRKC gene have been correlated with the development of medulloblastomas, secretory breast carcinomas, and other cancers. TRKC has several transcript variants that encode different isoforms (A= 94.4kDa, B= 68.4kDa, C= 92.8kDa, and D= 93.3kDa), and each isoform has unique signaling properties. Cell differentiation that plays a role in the development of proprioceptive neurons originates via signaling through TRKC. The TRKC protein contains two lg-like C2-type immunoglobulin-like domains, two leucine-rich repeats (LRR), one leucine rich repeat C-terminal (LRRCT) domain, and one protein kinase domain.

Synonyms: gp145(trkC); TRKC

Note: This TrkC antibody is useful for Immunohistochemistry-Paraffin and Immunocytochemistry /


Immunofluorescence. In ICC/IF, membrane staining was observed in neuro2a cells. In IHC-P, staining was observed in the membrane of mouse brain. Prior to immunostaining paraffin

tissues, antigen retrieval with sodium citrate buffer (pH 6.0) is recommended.

Protein Families: Druggable Genome, Protein Kinase, Transmembrane


Protein Pathways: Neurotrophin signaling pathway

Product images:



Immunocytochemistry/Immunofluorescence: TrkC Antibody TA336908 - TrkC antibody was tested in neuro2a cells with Dylight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and Dylight 550 (red).





Immunohistochemistry: TrkC Antibody TA336908 - TrkC antibody was tested in mouse brain using DAB with hematoxylin counterstain.