

Product datasheet for TA336940

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

EPCAM Mouse Monoclonal Antibody [Clone ID: 60N5D8]

Product data:

Product Type: Primary Antibodies

Clone Name: 60N5D8

Applications: FC, IHC, WB

Recommended Dilution: Immunohistochemistry: 1:10-1:500, Western Blot: 2-5 ug/ml, Immunohistochemistry-Paraffin:

1:200-1:400, Flow Cytometry: 0.2-1 ug/10^6 cells

Reactivity: Human Host: Mouse

Isotype: IgG2b, kappa
Clonality: Monoclonal

Immunogen: Within the range of amino acids 289-335 were used as the immunogen for the EpCAM

antibody.

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: lot specific

Purification: Protein G purified

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: epithelial cell adhesion molecule

Database Link: NP 002345

Entrez Gene 4072 Human

P16422

Synonyms: 4; DIAR5; EGP-2; EGP40; EGP314; ESA; HNPCC8; KS1; KSA; M4S1; MIC18; MK-1; TACSTD1;


TROP1

Protein Families: ES Cell Differentiation/IPS, Transmembrane






Product images:




Western Blot: EpCAM/TROP1 Antibody (60N5D8) TA336940 - Analysis of EpCAM in HT-29 lysate in the (A) absence and (B) presence of immunizing peptide using EpCAM antibody at 2 ug/ml. Goat anti-mouse Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this test.



Flow Cytometry: EpCAM/TROP1 Antibody (60N5D8) TA336940 - Analysis of HT-29 cells using EpCAM antibody at 0.2 ug/10^6 cells. this antibody was used for this test. Top picture shows no cell surface staining seen, bottom shows intracellular staining of EpCAM (red represents EpCAM antibody and green represents isotype control).





Immunohistochemistry-Paraffin: EpCAM/TROP1 Antibody (60N5D8) TA336940 - FFPE human colon stained with EpCAM antibody (1:200), peroxidase-conjugate and DAB chromogen.