

Product datasheet for **TA336953**

CILK1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1: 1000-1:5000
Reactivity:	Human, Bovine, Primate
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A portion of amino acid 380-430 of human ICK protein was used as the immunogen.
Formulation:	Store at -20C. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	intestinal cell (MAK-like) kinase
Database Link:	NP_057597 Entrez Gene 22858 Human Q9UPZ9
Synonyms:	ECO; LCK2; MRK
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Product images:



Western Blot: ICK Antibody TA336953 - Analysis of ICK in HeLa cell lysate. Lane 1 shows pre-immune sera. Lane 2 shows TA336953 tested at 1:5000 dilution.