

Product datasheet for **TA336963**

Rhox9 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western Blot: 1: 1000-1:2000
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Within the range of amino acids 82-112 of mouse Rhox-9 protein were used as the immunogen.
Formulation:	Store at -20C. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	reproductive homeobox 9
Database Link:	NP_076383 Entrez Gene 104384 Mouse
Synonyms:	1600026O01Rik; Gpbox; Psx2



[View online »](#)

Product images:



Western Blot: Rhox9 Antibody TA336963 - analysis of Rhox-9 in mouse placenta in the 1) absence and 2) presence of immunizing peptide using this antibody. I goat anti-rabbit Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this test.