

## Product datasheet for **TA336981**

### UVRAG Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western Blot: 0.5-1.0 ug/ml
Reactivity:	Human, Mouse, Rat, Canine, Chicken, Primate
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A portion of amino acids 375-425 of human UVRAG was used as the immunogen.
Formulation:	PBS, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	UV radiation resistance associated
Database Link:	<a href="#">NP_003360</a> <a href="#">Entrez Gene 78610 MouseEntrez Gene 308846 RatEntrez Gene 7405 Human Q9P2Y5</a>
Synonyms:	DHTX; p63; VPS38



[View online »](#)

## Product images:



Western Blot: UVRAG Antibody TA336981 - Analysis of UVRAG in human heart lysate in the 1) absence and 2) presence of immunizing peptide, 3) mouse heart and 4) rat heart using this antibody at 0 ug/ml, 0.5 ug/ml and 0.5 ug/ml, respectively.