

Product datasheet for **TA336984**

Syncollin (SYCN) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 0.05-2.0 ug/ml
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A portion of amino acids 50-100 of human syncollin was used as the immunogen.
Formulation:	PBS, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	syncollin
Database Link:	NP_001073937 Entrez Gene 342898 Human Q0VAF6
Synonyms:	INSSA1; SYL



[View online »](#)

Product images:



Western Blot: Syncollin Antibody TA336984 - Western blot analysis of syncollin in human pancreas lysate in the 1) absence and 2) presence of immunizing peptide and 3) rat pancreas using TA336984 at 0.05 ug/ml and 2 ug/ml, respectively.