

## **Product datasheet for TA337007**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Beta Arrestin 2 (ARRB2) Rabbit Polyclonal Antibody

**Product data:** 

**Product Type:** Primary Antibodies

**Applications:** IHC, Simple Western, WB

**Recommended Dilution:** Immunohistochemistry, Simple Western: 1:100, Western Blot: 0.25 - 0.5 ug/mL,

Immunohistochemistry-Paraffin: 5 ug/mL

**Reactivity:** Human, Mouse, Rat, Bovine, Canine, Equine, Primate

**Host:** Rabbit

Clonality: Polyclonal

**Immunogen:** A peptide sequence corresponding to an area between amino acids 1 and 50 of human

ARRB2 was used as the immunogen for the antibody, Gen Bank no. gb:ABG47460.1.

**Formulation:** PBS, 30% glycerol, 0.1% Sodium Azide. Store at -20C. Avoid freeze-thaw cycles.

**Concentration:** lot specific

**Purification:** Immunogen affinity purified

**Conjugation:** Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** arrestin beta 2 **Database Link:** NP 945355

Entrez Gene 25388 RatEntrez Gene 216869 MouseEntrez Gene 409 Human

P32121

Synonyms: ARB2; ARR2; BARR2


Protein Families: Druggable Genome

**Protein Pathways:** Chemokine signaling pathway, Endocytosis, MAPK signaling pathway, Olfactory transduction






## **Product images:**




Western Blot: beta-Arrestin 2 Antibody TA337007 - Analysis of ARRB2 in HUVEC cell lysate in the 1) absence and 2) presence of immunizing peptide using ARRB2 antibody at 0.25 ug/mL. Goat antirabbit Ig HRP secondary antibody and PicoTect ECL substrate solution were used.



Immunohistochemistry-Paraffin: beta-Arrestin 2 Antibody TA337007 - Human thyroid tissue stained with ARRB2 antibody at 5 ug/mL. Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM sodium citrate buffer, pH 6.0 for 10-20 min followed by cooling at RT for 20 min.





Simple Western: beta-Arrestin 2 Antibody TA337007 - Simple Western lane view shows a specific band for beta Arrestin in 0.5 mg/mL of Jurkat lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.