

Product datasheet for **TA337018**

RPS19BP1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western Blot: 2-4 ug/ml
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A portion of amino acids 1-50 of human AROS was used as the immunogen.
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ribosomal protein S19 binding protein 1
Database Link:	NP_919307 Entrez Gene 91582 Human Q86WX3
Synonyms:	AROS; S19BP



[View online »](#)

Product images:



Western Blot: RPS19BP1 Antibody TA337018 - Analysis of AROS in human liver lysate using this antibody at 3 ug/ml.