

Product datasheet for **TA337029**


ASAP3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	Immunohistochemistry, Immunohistochemistry-Paraffin: 5 ug/ml
Reactivity:	Human, Primate
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Within the range of amino acids 693-729 of human ASAP3 were used as the immunogen.
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ArfGAP with SH3 domain, ankyrin repeat and PH domain 3
Database Link:	NP_001137250 Entrez Gene 55616 Human Q8TDY4
Synonyms:	ACAP4; CENTB6; DDEFL1; UPLC1
Protein Pathways:	Endocytosis, Fc gamma R-mediated phagocytosis



[View online »](#)

Product images:

Immunohistochemistry-Paraffin: ACAP4 Antibody TA337029 - Analysis of ASAP3 in formalin-fixed, paraffin-embedded human brain tissue using an isotype control (top left) and this antibody (bottom left, right) at 5 ug/ml.