

Product datasheet for **TA337045**

BNC1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western Blot: 4-7 ug/ml, Knockdown Validated
Reactivity:	Human, Primate
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A portion of amino acids 330-380 of human Bnc 1 protein was used as the immunogen.
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Protein G purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	basonuclin 1
Database Link:	NP_001288135 Entrez Gene 646 Human Q01954
Synonyms:	BNC; BSN1; HsT19447



[View online »](#)

Product images:



Western Blot: BNC1 Antibody TA337045 - Analysis of BNC1 in human spleen lysate in the 1) absence and 2) presence of immunizing peptide, and 3) rat lysate using TA337045 at 7 ug/ml and 4 ug/ml, respectively.