

Product datasheet for **TA337082**


TLR6 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	Immunohistochemistry, Immunohistochemistry-Paraffin: 1: 100
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Within the range of amino acids 25-56 of human TLR6 protein were used as the immunogen.
Formulation:	Store at -20C. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	toll like receptor 6
Database Link:	NP_006059 Entrez Gene 10333 Human Q9Y2C9
Synonyms:	CD286
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Toll-like receptor signaling pathway



[View online »](#)

Product images:

Immunohistochemistry-Paraffin: TLR6 Antibody TA337082 - IHC analysis of TLR6 in paraffin-embedded human breast tissue using this antibody.