

Product datasheet for **TA337092**

FOXP3 Mouse Monoclonal Antibody [Clone ID: 16J4G6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	16J4G6
Applications:	WB
Recommended Dilution:	WB: 1-3 ug/ml
Reactivity:	Human
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Immunogen:	A synthetic peptide corresponding to the exon 2 deletion site was used to immunize mice.
Formulation:	PBS, 0.03% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Concentration:	lot specific
Purification:	Ammonium sulfate precipitation
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	forkhead box P3
Database Link:	NP_054728 Entrez Gene 50943 Human Q9BZS1
Synonyms:	AIID; DIETER; IPEX; JM2; PIDX; XPID
Protein Families:	Transcription Factors



[View online »](#)