

Product datasheet for **TA337339**

HEATR3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-HEATR3 antibody is: synthetic peptide directed towards the N-terminal region of Human HEATR3. Synthetic peptide located within the following region: KECSAGLDSNEMSLQEKKDQNRNSIENIANETVNLWNICECSSRAVSIF
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	35 kDa
Gene Name:	HEAT repeat containing 3
Database Link:	NP_891552 Entrez Gene 55027 Human Q7Z4Q2
Background:	The function of this protein remains unknown.
Synonyms:	SYO1
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Human: 100%; Rat: 92%; Rabbit: 92%; Bovine: 83%



[View online »](#)

Product images:



Host: Rabbit; Target Name: HEATR3; Sample Tissue: Fetal Kidney lysates; Antibody Dilution: 1.0 ug/ml