

Product datasheet for **TA337347**

HNRNPC Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-HNRNPC antibody is: synthetic peptide directed towards the middle region of Human HNRNPC. Synthetic peptide located within the following region: VDSLLENLEKIEKEQSKQAVEMKNDKSEEEQSSSVKKDETNVKMESEGG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	27 kDa
Gene Name:	heterogeneous nuclear ribonucleoprotein C (C1/C2)
Database Link:	NP_004491 Entrez Gene 15381 Mouse Entrez Gene 3183 Human P07910



[View online »](#)

Background:

This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene can act as a tetramer and is involved in the assembly of 40S hnRNP particles. Multiple transcript variants encoding at least two different isoforms have been described for this gene.

Synonyms:


C1; C2; HNRNP; HNRPC; SNRPC

Note:

Immunogen Sequence Homology: Dog: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Pig: 93%; Goat: 93%; Horse: 93%; Bovine: 93%; Rabbit: 93%; Guinea pig: 93%; Yeast: 85%

Protein Pathways:

Spliceosome

Product images:

Host: Rabbit; Target Name: HNRNPC; Sample Tissue: ACHN Whole Cell lysates; Antibody Dilution: 1.0 ug/ml There is BioGPS gene expression data showing that HNRNPC is expressed in ACHN