

Product datasheet for **TA337381**

KIAA0895L Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-KIAA0895L antibody is: synthetic peptide directed towards the N-terminal region of Human KIAA0895L. Synthetic peptide located within the following region: AGHIASKSPCMLVALRPTNMDRERDKFFQSHYTYPQFEYQEPMPPTAVLE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	44 kDa
Gene Name:	KIAA0895-like
Database Link:	NP_001035805 Entrez Gene 653319 Human Q68EN5
Background:	The function of this protein remains unknown.
Synonyms:	KIAA0895L
Note:	Immunogen Sequence Homology: Dog: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Pig: 86%; Mouse: 86%; Guinea pig: 86%



[View online »](#)

Product images:



Host: Rabbit; Target Name: KIAA0895L; Sample
Tissue: MCF7 whole cell lysates; Antibody
Dilution: 1.0 ug/ml