

Product datasheet for **TA337499**

OR4A15 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-OR4A15 antibody is: synthetic peptide directed towards the N-terminal region of Human OR4A15. Synthetic peptide located within the following region: KNNVTEFILLGLTQNPEGQKVLVFTLLIYMTIMGNLLIIVTIMASQSL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	38 kDa
Gene Name:	olfactory receptor family 4 subfamily A member 15
Database Link:	NP_001005275 Entrez Gene 81328 Human Q8NGL6



[View online »](#)

Background: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.


Synonyms: OR11-118

Note: Immunogen Sequence Homology: Rat: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Pig: 86%; Bovine: 86%; Dog: 85%; Horse: 85%; Guinea pig: 80%

Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

Product images:



Host: Rabbit; Target Name: OR4A15; Sample Tissue: HepG2 Whole Cell lysates; Antibody Dilution: 1.0 ug/ml