

## **Product datasheet for TA337681**

## **PIM1 Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

Applications: WB

Recommended Dilution: WB

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**Immunogen:** The immunogen for anti-PIM1 antibody: synthetic peptide directed towards the N terminal of

human PIM1. Synthetic peptide located within the following region: MLLSKINSLAHLRAAPCNDLHATKLAPGKEKEPLESQYQVGPLLGSGGFG

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Purification: Affinity Purified
Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 36 kDa

**Gene Name:** Pim-1 proto-oncogene, serine/threonine kinase

Database Link: NP 002639

Entrez Gene 18712 MouseEntrez Gene 5292 Human

P11309



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



**Background:** The protein encoded by this gene belongs to the Ser/Thr protein kinase family, and PIM

subfamily. This gene is expressed primarily in B-lymphoid and myeloid cell lines, and is overexpressed in hematopoietic malignancies and in prostate cancer. It plays a role in signal transduction in blood cells, contributing to both cell proliferation and survival, and thus provides a selective advantage in tumorigenesis. Both the human and orthologous mouse genes have been reported to encode two isoforms (with preferential cellular localization) resulting from the use of alternative in-frame translation initiation codons, the upstream non-

AUG (CUG) and downstream AUG codons.

Synonyms: PIM


**Note:** Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human:

100%; Mouse: 100%; Sheep: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase

**Protein Pathways:** Acute myeloid leukemia, Jak-STAT signaling pathway


## **Product images:**



WB Suggested Anti-PIM1 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control:

Human heart





Host: Rabbit; Target Name: PIM1; Sample Tissue: Human Fetal Lung; Antibody Dilution: 1.0ug/ml