


Product datasheet for **TA337712**

ZNF773 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF773 antibody: synthetic peptide directed towards the N terminal of human ZNF773. Synthetic peptide located within the following region: AAATLRDPAQQGYVTFEDVAVYFSQEEWRLLDDAQRLLYRNVMLENFTLL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	50 kDa
Gene Name:	zinc finger protein 773
Database Link:	NP_940944 Entrez Gene 374928 Human Q6PK81
Background:	May be involved in transcriptional regulation.
Synonyms:	ZNF419B
Note:	Immunogen Sequence Homology: Human: 100%; Pig: 79%
Protein Families:	Transcription Factors


[View online »](#)

Product images:

WB Suggested Anti-ZNF773 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control:
Human brain