

Product datasheet for **TA337773**

ZNF547 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF547 antibody: synthetic peptide directed towards the N terminal of human ZNF547. Synthetic peptide located within the following region: GTHPEQGLYTCPAHLHQHQKEQIREKLSRGDGGRRPTFVKNHRVHMAGKTF
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	38 kDa
Gene Name:	zinc finger protein 547
Database Link:	NP_775902 Entrez Gene 284306 Human Q8IVP9
Background:	May be involved in transcriptional regulation.
Synonyms:	FLJ31100
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZNF547 Antibody Titration:
2.5 ug/ml; Positive Control: Jurkat cell lysate