

Product datasheet for **TA337815**

C10orf46 (CACUL1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-CACUL1 antibody: synthetic peptide directed towards the middle region of human CACUL1. Synthetic peptide located within the following region: SELSEYAAQDQKFQRELIQNGFTRGDQSRKRAGDELAYNSSSACASSRGY
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	41 kDa
Gene Name:	CDK2 associated cullin domain 1
Database Link:	NP_722517 Entrez Gene 143384 Human Q86Y37
Background:	The specific function of the protein remains unknown.
Synonyms:	C10orf46; CAC1
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Guinea pig: 100%; Bovine: 93%



[View online »](#)

Product images:



WB Suggested Anti-C10orf46 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: HT1080 cell lysate