

Product datasheet for **TA337896**

ZNF396 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF396 antibody: synthetic peptide directed towards the middle region of human ZNF396. Synthetic peptide located within the following region: LQSLRPHEDEIKTTNVKSASRQKTSLGIELHCNVSNILHMNGSQSSTYRG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	38 kDa
Gene Name:	zinc finger protein 396
Database Link:	NP_665699 Entrez Gene 252884 Human Q96N95
Background:	ZNF396 acts as a DNA-dependent transcriptional repressor.
Synonyms:	ZSCAN14
Note:	Immunogen Sequence Homology: Human: 100%; Pig: 79%; Guinea pig: 79%; Rabbit: 77%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZNF396 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control:
Jurkat cell lysate