

## Product datasheet for **TA337981**

### PHF17 (JADE1) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PHF17 antibody: synthetic peptide directed towards the C terminal of human PHF17. Synthetic peptide located within the following region: EPFASLEQNREEAHRVSVRKQLQLEDEFYTFVNLDDVARALRLPEEVV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	58 kDa
Gene Name:	jade family PHD finger 1
Database Link:	<a href="#">NP_079176</a> <a href="#">Entrez Gene 79960 Human</a> <a href="#">Q6IE81</a>
Background:	PHF17 belongs to the JADE family. It contains 2 PHD-type zinc fingers. PHF17 is a transcriptional coactivator which seems to act by promoting acetylation of nucleosomal histone H4 by HTATIP. PHF17 promotes apoptosis. It may act as a renal tumor suppressor.
Synonyms:	PHF17
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%



[View online »](#)

Protein Families: Druggable Genome, Transcription Factors

### Product images:



WB Suggested Anti-PHF17 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control: Human brain