

## Product datasheet for **TA337997**


### **NDRG3 Rabbit Polyclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	WB
<b>Recommended Dilution:</b>	WB
<b>Reactivity:</b>	Human
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	The immunogen for anti-NDRG3 antibody is: synthetic peptide directed towards the C-terminal region of Human NDRG3. Synthetic peptide located within the following region: RLARSRTHSTSSSLGSGESPFERSVTSNQSDGTQESCESPDVLDLRHQTME
<b>Formulation:</b>	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
<b>Purification:</b>	Affinity Purified
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Predicted Protein Size:</b>	40 kDa
<b>Gene Name:</b>	NDRG family member 3
<b>Database Link:</b>	<a href="#">NP_071922</a> <a href="#">Entrez Gene 57446 Human</a> <a href="#">Q9UGV2</a>
<b>Background:</b>	The function of this protein remains unknown.
<b>Synonyms:</b>	FLJ13556
<b>Note:</b>	Immunogen Sequence Homology: Rat: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Dog: 93%; Pig: 93%; Guinea pig: 93%; Bovine: 92%; Horse: 86%



[View online »](#)

**Product images:**

WB Suggested Anti-NDRG3 Antibody; Titration: 1.0 ug/ml; Positive Control: OVCAR-3 Whole Cell  
NDRG3 is strongly supported by BioGPS gene expression data to be expressed in Human OVCAR3 cells