

Product datasheet for **TA338215**

UNC5C Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-UNC5C antibody: synthetic peptide directed towards the middle region of human UNC5C. Synthetic peptide located within the following region: VVVALFVYRKNHRDFESDIIDSSALNGGFQPVNIKAARQDLLAVPPDLTS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	103 kDa
Gene Name:	unc-5 netrin receptor C
Database Link:	NP_003719 Entrez Gene 8633 Human O95185



[View online »](#)

- Background:** This gene product belongs to the UNC-5 family of netrin receptors. Netrins are secreted proteins that direct axon extension and cell migration during neural development. They are bifunctional proteins that act as attractants for some cell types and as repellents for others, and these opposite actions are thought to be mediated by two classes of receptors. The UNC-5 family of receptors mediate the repellent response to netrin; they are transmembrane proteins containing 2 immunoglobulin (Ig)-like domains and 2 type I thrombospondin motifs in the extracellular region. [provided by RefSeq, Jul 2008]
- Synonyms:** UNC5H3
- Note:** Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Zebrafish: 86%
- Protein Families:** Druggable Genome, Transmembrane
- Protein Pathways:** Axon guidance
- Product images:**



WB Suggested Anti-UNC5C Antibody Titration:
0.2-1 ug/ml; Positive Control: COLO205 cell lysate