

Product datasheet for **TA338646**

PVRIG Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PVRIG antibody: synthetic peptide directed towards the middle region of human PVRIG. Synthetic peptide located within the following region: ACGSLPPSSDPGLSAPPTPAPILRADLAGILGVSGVLLFGCVYLLHLLRR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34 kDa
Gene Name:	poliovirus receptor related immunoglobulin domain containing
Database Link:	NP_076975 Entrez Gene 79037 Human Q6DKI7
Background:	The exact function of PVRIG remains unknown.
Synonyms:	C7orf15
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transmembrane


[View online »](#)

Product images:



WB Suggested Anti-PVRIG Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive Control: Human brain