

Product datasheet for **TA338666**

Ogfod3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-RGD1305007 antibody is: synthetic peptide directed towards the N-terminal region of Rat RGD1305007. Synthetic peptide located within the following region: RFE _G CTPRKCGRGVTDIVITREEAEQIRRIAEGKLSLGGSDGGASILD _{LH}
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34 kDa
Gene Name:	2-oxoglutarate and iron-dependent oxygenase domain containing 3
Database Link:	NP_001013998 Entrez Gene 303749 Rat Q5M843
Background:	The function of this protein remains unknown.
Synonyms:	C17orf101
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Guinea pig: 100%; Rabbit: 93%; Zebrafish: 93%



[View online »](#)

Product images:



Host: Rabbit; Target Name: RGD1305007; Sample Tissue: Rat Testis lysates; Antibody Dilution: 1.0 ug/ml