

Product datasheet for TA338877


C11orf65 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-C11orf65 antibody: synthetic peptide directed towards the N terminal of human C11orf65. Synthetic peptide located within the following region: MPWKEESEFTKQDKAARVIQAWKSFLNVAIFQHFSLIDLRQGEPRQI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	37 kDa
Gene Name:	chromosome 11 open reading frame 65
Database Link:	NP_689800 Entrez Gene 160140 Human Q8NCR3
Background:	The specific function of C11orf65 is not yet known.
Synonyms:	MGC33948
Note:	Immunogen Sequence Homology: Rat: 100%; Human: 100%; Mouse: 100%; Pig: 93%; Guinea pig: 93%; Horse: 86%; Rabbit: 86%


[View online »](#)

Product images:



WB Suggested Anti-C11orf65 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: MCF7 cell lysate