

Product datasheet for **TA339125**


ZNF230 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF230 antibody: synthetic peptide directed towards the N terminal of human ZNF230. Synthetic peptide located within the following region: IHTGQKPSQNGKCKQSFSDVAIFDPPQQFHSGEKSHTCNECGKSFYISA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54 kDa
Gene Name:	zinc finger protein 230
Database Link:	NP_006291 Entrez Gene 7773 Human Q9UIE0
Background:	May be involved in transcriptional regulation.
Synonyms:	FDZF2
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

Product images:

WB Suggested Anti-ZNF230 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
Control: Jurkat cell lysate