

Product datasheet for **TA339164**

Rps12 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-Rps12 antibody is: synthetic peptide directed towards the C-terminal region of Rps12. Synthetic peptide located within the following region: NKKLG EWVGLCKIDREGKPRKVVGCSCVVKDYGKESQAKDVIEEYFKCK
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	14 kDa
Gene Name:	ribosomal protein S12
Database Link:	NP_035425 Entrez Gene 6206 Human Entrez Gene 20042 Mouse
Background:	The function of this protein remains unknown.
Synonyms:	OTTHUMP00000017230
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Sheep: 100%; Bovine: 100%; Rabbit: 100%; Zebrafish: 100%; Guinea pig: 100%



[View online »](#)

Product images:



WB Suggested Anti-Rps12 Antibody; Titration: 1.0 ug/ml; Positive Control: Mouse Pancreas