

Product datasheet for **TA339228**

ARID5A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ARID5A antibody: synthetic peptide directed towards the N terminal of human ARID5A. Synthetic peptide located within the following region: MAAPVKGNRKQSTEGDALDPPASPKPAGKQNGIQNPISLEDSPEAGGERE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	64 kDa
Gene Name:	AT-rich interaction domain 5A
Database Link:	NP_997646 Entrez Gene 10865 Human Q03989
Background:	Members of the ARID protein family, including ARID5A, have diverse functions but all appear to play important roles in development, tissue-specific gene expression, and regulation of cell growth. Members of the ARID protein family, including ARID5A, have diverse functions but all appear to play important roles in development, tissue-specific gene expression, and regulation of cell growth (Patsialou et al., 2005 [PubMed 15640446]). [supplied by OMIM]




[View online »](#)

Synonyms: MRF-1; MRF1; RP11-363D14

Note: Immunogen Sequence Homology: Human: 100%

Product images:



WB Suggested Anti-ARID5A Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: COLO205 cell lysate