

Product datasheet for **TA339445**

ZNF195 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF195 antibody: synthetic peptide directed towards the middle region of human ZNF195. Synthetic peptide located within the following region: QCEECGKVFRTCSSLSNHKRTTHSEEKPYTCEEKGNIFKQLSDLTKHKKTH
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	72 kDa
Gene Name:	zinc finger protein 195
Database Link:	NP_001123992 Entrez Gene 7748 Human O14628



[View online »](#)


Background: This gene encodes a protein belonging to the Krueppel C2H2-type zinc-finger protein family. These family members are transcription factors that are implicated in a variety of cellular processes. This gene is located near the centromeric border of chromosome 11p15.5, next to an imprinted domain that is associated with maternal-specific loss of heterozygosity in Wilms' tumors. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]

Synonyms: HRF1; ZNFP104

Note: Immunogen Sequence Homology: Human: 100%

Protein Families: Transcription Factors

Product images:



WB Suggested Anti-ZNF195 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:1562500; Positive
Control: THP-1 cell lysate