

## Product datasheet for **TA339461**

### ZNF225 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF225 antibody: synthetic peptide directed towards the N terminal of human ZNF225. Synthetic peptide located within the following region: QDDMPCQVDAGLSIIHVKTETSEGRTCKKSFSDVSVLDLHQQLSREKSH
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	82 kDa
Gene Name:	zinc finger protein 225
Database Link:	<a href="#">NP_037494</a> <a href="#">Entrez Gene 7768 Human</a> <a href="#">Q9UK10</a>
Background:	May be involved in transcriptional regulation.
Synonyms:	MGC119735
Note:	Immunogen Sequence Homology:
Protein Families:	Transcription Factors



[View online »](#)

## Product images:



WB Suggested Anti-ZNF225 Antibody Titration:  
0.3125 ug/ml; ELISA Titer: 1:312500; Positive  
Control: HepG2 cell lysate