

Product datasheet for **TA339475**


ZNF510 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF510 antibody: synthetic peptide directed towards the N terminal of human ZNF510. Synthetic peptide located within the following region: SEEEFSNQSHPKDYRGDDLIQNKIKDKHLEQAICINNKTLTTEEEKVL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	79 kDa
Gene Name:	zinc finger protein 510
Database Link:	NP_055745 Entrez Gene 22869 Human Q9Y2H8
Background:	May be involved in transcriptional regulation.
Synonyms:	KIAA0972; MGC33740
Note:	Immunogen Sequence Homology: Human: 100%; Pig: 75%
Protein Families:	Transcription Factors



[View online »](#)

Product images:

WB Suggested Anti-ZNF510 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
Control: Human Thymus