

Product datasheet for **TA339519**

ZNF490 Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Recommended Dilution: | WB |
| Reactivity: | Human |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen for anti-ZNF490 antibody: synthetic peptide directed towards the N terminal of human ZNF490. Synthetic peptide located within the following region: DEHKNQGRNLRSPMVEALCENKEDCPCGKSTSQIPDLNTNLETPTGLKPC |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i> |
| Concentration: | lot specific |
| Purification: | Protein A purified |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Predicted Protein Size: | 61 kDa |
| Gene Name: | zinc finger protein 490 |
| Database Link: | NP_065765 Entrez Gene 57474 Human Q9ULM2 |
| Background: | May be involved in transcriptional regulation. |
| Synonyms: | KIAA1198 |
| Note: | Immunogen Sequence Homology: Human: 100% |
| Protein Families: | Transcription Factors |



[View online »](#)

Product images:



WB Suggested Anti-ZNF490 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
Control: OVCAR-3 cell lysate