

Product datasheet for **TA339576**

FCRL4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-FCRL4 antibody: synthetic peptide directed towards the N terminal of human FCRL4. Synthetic peptide located within the following region: FKGERVTLTCNGFQFYATEKTTWYHRHYWGEKLTLPGNTLEVRESGLYR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57 kDa
Gene Name:	Fc receptor like 4
Database Link:	NP_112572 Entrez Gene 83417 Human Q96PJ5
Background:	FCRL4 may function as an inhibitor of the B cell receptor signaling and in the B cell-mediated immune response.
Synonyms:	CD307d; FCRH4; IGFP2; IRTA1
Note:	Immunogen Sequence Homology: Human: 100%; Pig: 85%
Protein Families:	Transmembrane



[View online »](#)

Product images:

WB Suggested Anti-FCRL4 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 2500; Positive Control: 293T cell lysate