

Product datasheet for **TA339591**

FAM62C (ESYT3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-FAM62C antibody: synthetic peptide directed towards the middle region of human FAM62C. Synthetic peptide located within the following region: EVFEFMVYEVPGQDLEVDLYDEDTDRDDFLGSLQICLGDMVTNRVDEWF
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	100 kDa
Gene Name:	extended synaptotagmin protein 3
Database Link:	NP_114119 Entrez Gene 83850 Human A0FGR9
Background:	ESYT3 belongs to the extended synaptotagmin family. It is a single-pass membrane protein, and contains 3 C2 domains. ESYT3 may play a role as calcium-regulated intrinsic membrane protein.
Synonyms:	CHR3SYT; E-Syt3; FAM62C




[View online »](#)

Note: Immunogen Sequence Homology: Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Guinea pig: 100%; Dog: 93%; Rabbit: 93%; Horse: 92%; Bovine: 92%; Zebrafish: 77%

Protein Families: Transmembrane

Product images:



WB Suggested Anti-FAM62C Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: Human Stomach