

Product datasheet for **TA339646**

ZMYND17 (MSS51) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZMYND17 antibody: synthetic peptide directed towards the C terminal of human ZMYND17. Synthetic peptide located within the following region: SAHRGLYHDFWEEQVETGQTHHPDLVAAFHPGFHSSPDLMEAWLPTLLLL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	51 kDa
Gene Name:	MSS51 mitochondrial translational activator
Database Link:	NP_001019764 Entrez Gene 118490 Human Q4VC12
Background:	ZMYND17 contains 1 MYND-type zinc finger. The exact function of ZMYND17 remains unknown.
Synonyms:	ZMYND17
Note:	Immunogen Sequence Homology: Human: 100%; Dog: 86%; Pig: 79%; Horse: 79%; Bovine: 79%; Rabbit: 79%; Guinea pig: 79%



[View online »](#)

Product images:



WB Suggested Anti-MSS51 Antibody Titration:
0.2-1 ug/ml; Positive Control: HepG2 cell lysate