

Product datasheet for **TA339739**

NRIP3 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-NRIP3 antibody: synthetic peptide directed towards the middle region of human NRIP3. Synthetic peptide located within the following region: ALVDTGCLYNLISLACVDRLGLKEHVKSHKHEGEKLSLPRHLKVVQIEH
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	27 kDa
Gene Name:	nuclear receptor interacting protein 3
Database Link:	NP_065696 Entrez Gene 56675 Human Q9NQ35
Background:	The exact functions of NRIP3 remain unknown.
Synonyms:	C11orf14; NY-SAR-105
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Guinea pig: 100%; Bovine: 92%



[View online »](#)

Protein Families: Druggable Genome

Product images:



WB Suggested Anti-NRIP3 Antibody Titration: 0.2-1 ug/ml; Positive Control: 721_B cell lysate