

Product datasheet for **TA339774**


ZMAT1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZMAT1 antibody: synthetic peptide directed towards the middle region of human ZMAT1. Synthetic peptide located within the following region: RMCEQRFSEASQTYQRPYHISPVESQLPQWLPTHSKRTYDSFQDELEDY
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55 kDa
Database Link:	NP_115817 Entrez Gene 84460 Human
Background:	The function of ZMAT1 has not been determined.
Synonyms:	KIAA1789; matrin type 1; MGC176597; MGC176728; OTTHUMP00000023703; zinc finger
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Horse: 100%; Human: 100%; Rat: 92%; Mouse: 92%; Bovine: 92%; Rabbit: 92%; Guinea pig: 92%



[View online »](#)

Product images:

WB Suggested Anti-ZMAT1 Antibody Titration:
0.2-1 ug/ml; Positive Control: HepG2 cell lysate