


Product datasheet for **TA339786**

ZNF618 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF618 antibody: synthetic peptide directed towards the N terminal of human ZNF618. Synthetic peptide located within the following region: TLGSYECGICGKKYKYNCFCQTHVRAHRDTEATSGEGASQSNNFRYTCDI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	41 kDa
Gene Name:	zinc finger protein 618
Database Link:	NP_588615 Entrez Gene 114991 Human Q5T7W0
Background:	ZNF618 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 4 C2H2-type zinc fingers. ZNF618 may be involved in transcriptional regulation.
Synonyms:	FP13169; NEDD10
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Mouse: 93%


[View online »](#)

Product images:


WB Suggested Anti-ZNF618 Antibody Titration:
0.2-1 ug/ml; Positive Control: NTERA2 cell lysate