

Product datasheet for **TA340208**


CTRP4 (C1QTNF4) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-C1QTNF4 antibody: synthetic peptide directed towards the middle region of human C1QTNF4. Synthetic peptide located within the following region: RRGDAVWLLSHDHDGYGAYSNHGKYITFSGFLVYPDLAPAAPPGLGASEL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	35 kDa
Gene Name:	C1q and tumor necrosis factor related protein 4
Database Link:	NP_114115 Entrez Gene 114900 Human Q9BXJ3
Background:	The specific function of C1QTNF4 is not yet known.
Synonyms:	CTRP4; ZACRP4
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Guinea pig: 100%; Zebrafish: 83%
Protein Families:	Secreted Protein



[View online »](#)

Product images:

WB Suggested Anti-C1QTNF4 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 12500; Positive
Control: Transfected 293T