

Product datasheet for **TA340209**

CTRP7 (C1QTNF7) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-C1QTNF7 antibody: synthetic peptide directed towards the middle region of human C1QTNF7. Synthetic peptide located within the following region: SIVLKSAFSVGIITTSYPEERLPIIFNKVLFNEGEHYNPATGKFICAFPGI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	31 kDa
Gene Name:	C1q and tumor necrosis factor related protein 7
Database Link:	NP_114117 Entrez Gene 114905 Human Q9BXJ2
Background:	The specific function of C1QTNF7 is not yet known.
Synonyms:	CTRP7; ZACRP7
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 93%; Zebrafish: 86%



[View online »](#)

Protein Families: Secreted Protein

Product images:



WB Suggested Anti-C1QTNF7 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: COLO205 cell lysate