

Product datasheet for **TA340222**

LRRC42 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-LRRC42 antibody: synthetic peptide directed towards the N terminal of human LRRC42. Synthetic peptide located within the following region: LIGFPEQIAEKLFSAAEARQKFTEPGAGLRALQKFTEAYGSLVLCSLCLR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	48 kDa
Gene Name:	leucine rich repeat containing 42
Database Link:	NP_443172 Entrez Gene 115353 Human Q9Y546
Background:	The specific function of this protein remains unknown.
Synonyms:	dj167A19.4
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 93%



[View online »](#)

Product images:

WB Suggested Anti-LRRC42 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: HT1080 cell lysate LRRC42 is supported
by BioGPS gene expression data to be expressed
in HT1080