

## **Product datasheet for TA340639**

## **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **GPR124 (ADGRA2) Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

Applications: IHC

**Recommended Dilution:** IHC-P (4 µg/ml)

**Reactivity:** Human (Predicted: Mouse, Rat, Hamster, Chicken, Xenopus)

**Host:** Rabbit

Clonality: Polyclonal

Immunogen: GPR124 / TEM5 antibody was raised against synthetic 17 amino acid peptide from N-terminal

extracellular domain of human GPR124. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Marmoset, Dog, Bat, Elephant, Panda, Horse, Rabbit, Pig, Opossum (100%); Mouse, Rat, Hamster, Turkey, Chicken, Xenopus (94%); Monkey, Bovine,

Platypus, Lizard (82%).

**Formulation:** PBS, 0.1% sodium azide.

**Concentration:** lot specific

Purification: Immunoaffinity purified

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: adhesion G protein-coupled receptor A2

Database Link: NP 116166

Entrez Gene 78560 MouseEntrez Gene 100363275 RatEntrez Gene 25960 Human

Q96PE1


**Synonyms:** GPR124; TEM5

**Protein Families:** Druggable Genome, Transmembrane






## **Product images:**




Anti-GPR124 / TEM5 antibody IHC of human Prostate, Carcinoma Endothelial Cells. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.




Anti-GPR124 antibody IHC of human colon carcinoma, neoplastic cells. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.





Anti-GPR124 / TEM5 antibody IHC of human Lung, Adenocarcinoma. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.



Anti-GPR124 / TEM5 antibody IHC of human Skin, Kaposis Sarcoma capillary. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.