

Product datasheet for **TA340890**


MRGPRX4 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (6 - 28 µg/ml)
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	MRGPRX4 / MRGX4 antibody was raised against synthetic 12 amino acid peptide from N-terminal extracellular domain of human MRGPRX4 / MRGX4. Percent identity with other species by BLAST analysis: Human (100%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	MAS related GPR family member X4
Database Link:	NP_473373 Entrez Gene 117196 Human Q96LA9
Synonyms:	GPCR; MRGX4; SNSR6
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)

Product images:

Anti-MRGPRX4 / MRGX4 antibody IHC of human Lung, Non-Small Cell Carcinoma.
Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.



Anti-MRGPRX4 / MRGX4 antibody IHC of human spinal cord, dorsal root ganglion.
Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.