

# Product datasheet for TA341067

## NUR77 (NR4A1) Rabbit Polyclonal Antibody

### **Product data:**

#### OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (16 μg/ml)
Reactivity:	Human (Predicted: Monkey, Bovine, Dog, Horse)
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	NR4A1 / NUR77 antibody was raised against synthetic 15 amino acid peptide from N-terminal to DNA binding domain of human NUR77. Percent identity with other species by BLAST analysis: Human, Gorilla, Orangutan, Gibbon, Monkey, Mouse, Rat, Hamster (100%); Marmoset, Elephant, Panda, Bovine, Dog, Horse (93%); Bat, Rabbit (87%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	nuclear receptor subfamily 4 group A member 1
Database Link:	<u>NP_002126</u> <u>Entrez Gene 403897 DogEntrez Gene 696797 MonkeyEntrez Gene 3164 Human</u> <u>P22736</u>
Synonyms:	GFRP1; HMR; N10; NAK-1; NGFIB; NP10; NUR77; TR3
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	MAPK signaling pathway



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US 

### **Product images:**



Anti-NUR77 antibody IHC of human ovary. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US