

Product datasheet for **TA341108**


UCHL3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (5 - 10 µg/ml)
Reactivity:	Human (Predicted: Monkey, Bovine, Horse, Pig, Rabbit)
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	UCHL3 antibody was raised against synthetic 18 amino acid peptide from internal region of human UCHL3. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Mouse, Rat, Dog, Elephant, Opossum, Turkey, Chicken, Lizard (100%); Orangutan, Monkey, Bovine, Horse, Rabbit, Pig, Salmon (94%); Pufferfish, Zebrafish, Stickleback, Water flea (89%); Xenopus (83%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ubiquitin C-terminal hydrolase L3
Database Link:	NP_005993 Entrez Gene 697276 Monkey Entrez Gene 7347 Human P15374
Synonyms:	UCH-L3
Protein Families:	Druggable Genome, Protease



[View online »](#)

Product images:

Human Skeletal Muscle: Formalin-Fixed, Paraffin-Embedded (FFPE)