

Product datasheet for **TA341150**


Fibroblast activation protein, alpha (FAP) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (5 µg/ml)
Reactivity:	Human (Predicted: Bovine, Pig)
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	FAP Alpha antibody was raised against synthetic 17 amino acid peptide from extracellular domain of human FAP. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Horse (100%); Panda, Bovine, Pig (94%); Elephant, Dog (88%); Rat, Bat, Platypus (82%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	fibroblast activation protein alpha
Database Link:	NP_004451 Entrez Gene 2191 Human Q12884
Synonyms:	DPPIV; FAPA; SIMP
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)

Product images:

Anti-FAP antibody IHC of human skin, hair follicle. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.