

Product datasheet for **TA341191**


ADAMTS4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (8 µg/ml)
Reactivity:	Human (Predicted: Monkey)
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	ADAMTS4 antibody was raised against synthetic 19 amino acid peptide from internal region of human ADAMTS4. Percent identity with other species by BLAST analysis: Human, Orangutan (100%); Gibbon, Marmoset (95%); Monkey, Elephant, Horse (89%); Bovine, Dog, Panda, Opossum (84%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ADAM metallopeptidase with thrombospondin type 1 motif 4
Database Link:	NP_005090 Entrez Gene 719927 Monkey Entrez Gene 9507 Human O75173
Synonyms:	ADAMTS-2; ADAMTS-4; ADMP-1
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane



[View online »](#)

Product images:

Anti-ADAMTS4 antibody IHC of human heart. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.